

Indications of EyeBrid lenses



Indications of **all** EyeBrid Silicone lenses

Wearers of RGPs

Problems relating to comfort/centration/stability

Occasional wear

Keratoconus and other irregular corneas

including highly complex irregular corneas

Astigmatisms

that can potentially be fitted with RGPs

High ametropias

And Piggyback refits, 1st generation hybrid lens refits, Unilateral wear

Particular indications of EyeBrid Silicone **Multifocal**

EyeBrid Silicone **Multifocal** lenses are indicated for **presbyopia** with the indications as listed above, and for regular and irregular **astigmatism** in particular.

As a **general rule**

EyeBrid Silicone and EyeBrid Silicone Multifocal lenses are:

- **Indicated** for RGP refits and potential RGP wearers
- **Not recommended** for former wearers of soft lenses

Any questions? Contact technicalsupport@eyebbridsilicone.com